

CAPG Human

Capping Protein Gelsolin-Like Human Recombinant

RCP0333

Product Overview

Name	CAPG Human
Accession(Primary)	P40121
Immunogen	Anti-human CAPNS1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CAPNS1 protein 84-268 amino acids purified from E. coli.
Description	Capping Protein Gelsolin-Like Human Recombinant

Target information (P40121)

Gene Name	
About Gene Symbol	CAPG (Capping Actin Protein, Gelsolin Like) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.