

CBX5 Human

Chromobox Homolog 5 Human Recombinant

RCP0347

Product Overview

Name	CBX5 Human
Accession(Primary)	P45973
Description	Chromobox Homolog 5 Human Recombinant

Target information (P45973)

Gene Name

About Gene Symbol

CBX5 (Chromobox) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics and adaptive immunity. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.