

CHD4 Human

Chromodomain Helicase DNA Binding Protein 4 Human Recombinant

RCP0419

Product Overview

Name	CHD4 Human
Accession(Primary)	Q14839
Immunogen	Anti-human CHGA mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human CHGA amino acids 19-131 purified from E. coli.
Description	Chromodomain Helicase DNA Binding Protein 4 Human Recombinant

Target information (Q14839)

Gene Name	
About Gene Symbol	CHD4 (Chromodomain Helicase DNA Binding Protein 4) is an epigenetic and nuclear regulatory protein involved in adaptive immunity. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.