

CHGB Human

Chromogranin B Human Recombinant

RCP0424

Product Overview

Name	CHGB Human
Accession(Primary)	P05060
Description	Chromogranin B Human Recombinant

Target information (P05060)

Gene Name

About Gene Symbol

References for protein:

Title: Secretoneurin promotes neuroprotection and neuronal plasticity via the Jak2/Stat3 pathway in murine models of stroke. Publication: The Journal of Clinical Investigation, Volume 118, Number 1, January 2008. Link: CHGB prospec publication Applications: CHGB acts directly on neurons after hypoxia and ischemic insult to further their survival by activating the Jak2/Stat3 pathway.