

CIRBP Human

Cold Inducible RNA Binding Protein Human Recombinant

RCP0440

Product Overview

Name	CIRBP Human
Accession(Primary)	Q14011
Immunogen	Anti-human CISD1 mAb, clone PAT1A8A, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CISD1 protein 32-108 amino acids purified from E. coli.
Description	Cold Inducible RNA Binding Protein Human Recombinant

Target information (Q14011)

Gene Name

About Gene Symbol

CIRBP (Cold Inducible RNA Binding Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.