

CLIC4 Human

Chloride Intracellular Channel 4 Human Recombinant

RCP0459

Product Overview

Name	CLIC4 Human
Accession(Primary)	Q9Y696
Description	Chloride Intracellular Channel 4 Human Recombinant

Target information (Q9Y696)

Gene Name

About Gene Symbol

CLIC4 (Chloride Intracellular Channel 4) is a membrane transporter involved in membrane transport. It is an important research target in pharmacology, drug resistance, and membrane biology and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.