

DDX56 Human

DEAD Box Protein 56 Human Recombinant

RCP0631

Product Overview

Name	DDX56 Human
Accession(Primary)	Q9NY93
Immunogen	Anti-human DECR1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human DECR1 protein 35-335 amino acids purified from E. coli.
Description	DEAD Box Protein 56 Human Recombinant

Target information (Q9NY93)

Gene Name

About Gene Symbol