

DKK2 Human

Dickkopf-Related Protein 2 Human Recombinant

RCP0644

Product Overview

Name	DKK2 Human
Accession(Primary)	Q9UBU2
Immunogen	Anti-human DKK3 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human DKK3 protein (22-350 amino acids purified from E. coli.
Description	Dickkopf-Related Protein 2 Human Recombinant

Target information (Q9UBU2)

Gene Name

About Gene Symbol