

Dnak SBD C-terminus E.Coli

Dnak Substrate Binding Domain C-terminal E.Coli Recombinant

RCP0670

Product Overview

Name	Dnak SBD C-terminus E.Coli
Immunogen	Anti-human DNAL1 mAb, clone PAT29E4AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human DNAL1 protein 1-190 amino acids purified from E. coli.
Description	Dnak Substrate Binding Domain C-terminal E.Coli Recombinant