

DYNLL1 Human

Dynein Light Chain LC8 Type-1 Human Recombinant

RCP0689

Product Overview

Name	DYNLL1 Human
Accession(Primary)	P63167
Description	Dynein Light Chain LC8 Type-1 Human Recombinant

Target information (P63167)

Gene Name

About Gene Symbol

DYNLL1 (Dynein Light Chain LC8-Type 1) is a biosynthetic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.