

EIF4E Human

Eukaryotic Translation Initiation Factor 4E Human Recombinant

RCP0741

Product Overview

Name	EIF4E Human
Accession(Primary)	P06730
Description	Eukaryotic Translation Initiation Factor 4E Human Recombinant

Target information (P06730)

Gene Name

About Gene Symbol

EIF4E (Eukaryotic Translation Initiation Factor 4E) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.