

EIF4H Human

Eukaryotic Translation Initiation Factor 4H Human Recombinant

RCP0746

Product Overview

Name	EIF4H Human
Accession(Primary)	Q15056
Immunogen	Anti-human EIF5A mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human EIF5A amino acids 1-154 purified from E. coli.
Description	Eukaryotic Translation Initiation Factor 4H Human Recombinant

Target information (Q15056)

Gene Name

About Gene Symbol