

ELP4 Human

Elongator Acetyltransferase Complex Subunit 4 Human Recombinant

RCP0752

Product Overview

Name	ELP4 Human
Accession(Primary)	Q96EB1
Description	Elongator Acetyltransferase Complex Subunit 4 Human Recombinant

Target information (Q96EB1)

Gene Name

About Gene Symbol

ELP4 (Elongator Acetyltransferase Complex Subunit 4) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.