

EPHA2 Human, sf9

EPH Receptor A2 Human Recombinant, sf9

RCP0763

Product Overview

Name	EPHA2 Human, sf9
Accession(Primary)	P29317
Immunogen	Anti-human EPHA2 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human EPHA2 protein 27-537 amino acids purified from E. coli.
Description	EPH Receptor A2 Human Recombinant, sf9

Target information (P29317)

Gene Name**About Gene Symbol**