

FABP2 Human, His

Fatty Acid Binding Protein 2 Human Recombinant, His Tag

RCP0809

Product Overview

Name	FABP2 Human, His
Accession(Primary)	P12104
Immunogen	Anti-human FABP3 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human FABP3 protein 1-133 amino acids purified from E. coli.
Description	Fatty Acid Binding Protein 2 Human Recombinant, His Tag

Target information (P12104)

Gene Name**About Gene Symbol**