

FIMC E.Coli

Chaperone Protein fimC E.Coli Recombinant

RCP0857

Product Overview

Name	FIMC E.Coli
Accession(Primary)	P31697
Immunogen	Anti-human FIS1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human FIS1 protein 1-122 amino acids purified from E. coli.
Description	Chaperone Protein fimC E.Coli Recombinant

Target information (P31697)

Gene Name

About Gene Symbol