

FSBP Human

Fibrinogen Silencer Binding Protein Human Recombinant

RCP0876

Product Overview

Name	FSBP Human
Accession(Primary)	O95073
Immunogen	Anti-human FSCN1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human FSCN1 protein 1-493 amino acids purified from E. coli.
Description	Fibrinogen Silencer Binding Protein Human Recombinant

Target information (O95073)

Gene Name

About Gene Symbol