

GAP43 Human

Growth Associated Protein 43 Human Recombinant

RCP0898

Product Overview

Name	GAP43 Human
Accession(Primary)	P17677
Immunogen	Anti-human GAPDH mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human GAPDH amino acids 1-335 purified from E. coli.
Description	Growth Associated Protein 43 Human Recombinant

Target information (P17677)

Gene Name

About Gene Symbol

GAP43 (Growth Associated Protein 43) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.