

GNB2L1 Human

Guanine Nucleotide Binding Protein beta polypeptide 2-like 1 Human Recombinant

RCP0940

Product Overview

Name	GNB2L1 Human
Accession(Primary)	P63244
Description	Guanine Nucleotide Binding Protein beta polypeptide 2-like 1 Human Recombinant

Target information (P63244)

Gene Name

About Gene Symbol

RACK1 (Receptor For Activated C Kinase 1) is a membrane receptor involved in translational regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.