

GPC4 511 aa Human

Glypican-4 511 aa Human Recombinant

RCP0957

Product Overview

Name	GPC4 511 aa Human
Accession(Primary)	O75487
Immunogen	Anti-human GPC4 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human GPC4 protein 401-529 amino acids purified from E. coli.
Description	Glypican-4 511 aa Human Recombinant

Target information (O75487)

Gene Name**About Gene Symbol**