

HINT1 Human

Histidine Triad Nucleotide Binding Protein 1 Human Recombinant

RCP1019

Product Overview

Name	HINT1 Human
Accession(Primary)	P49773
Immunogen	Anti-human HINT2 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human HINT2 protein 18-163 amino acids purified from E. coli.
Description	Histidine Triad Nucleotide Binding Protein 1 Human Recombinant

Target information (P49773)

Gene Name

About Gene Symbol