

HMGA2 Human

High Mobility Group AT-Hook 2 Human Recombinant

RCP1032

Product Overview

Name	HMGA2 Human
Accession(Primary)	P52926
Immunogen	Anti-human HMGB1 mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human HMGB1 amino acids 1-215 purified from High Five (Trichoplusia ni) insect cells.
Description	High Mobility Group AT-Hook 2 Human Recombinant

Target information (P52926)

Gene Name

About Gene Symbol