

HMGN3 Human

High Mobility Group Nucleosomal Binding Domain 3 Human Recombinant

RCP1038

Product Overview

Name	HMGN3 Human
Accession(Primary)	Q15651
Immunogen	Anti-human HMOX1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human HMOX1 amino acids 1-266 purified from E. coli.
Description	High Mobility Group Nucleosomal Binding Domain 3 Human Recombinant

Target information (Q15651)

Gene Name

About Gene Symbol