

HNRNPK Human

Heterogeneous Nuclear Ribonucleoprotein K Human Recombinant

RCP1044

Product Overview

Name	HNRNPK Human
Accession(Primary)	P61978
Immunogen	Anti-human HOMER1 mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human HOMER1 amino acids 1-354 purified from E. coli.
Description	Heterogeneous Nuclear Ribonucleoprotein K Human Recombinant

Target information (P61978)

Gene Name	
About Gene Symbol	HNRNPK (Heterogeneous Nuclear Ribonucleoprotein K) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.