

ICOSLG Human

Inducible T-Cell Costimulator Ligand Human Recombinant

RCP1105

Product Overview

Name	ICOSLG Human
Accession(Primary)	O75144
Immunogen	Anti-human ICT1 mAb, clone PAT1E9A, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ICT1 protein 30-206 amino acids purified from E. coli.
Description	Inducible T-Cell Costimulator Ligand Human Recombinant

Target information (O75144)

Gene Name

About Gene Symbol