

ING1 Human

Inhibitor of Growth Family, Member 1 Human Recombinant

RCP1127

Product Overview

Name	ING1 Human
Accession(Primary)	Q9UK53
Immunogen	Anti-human ING2 mAb, clone PAT39E5AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ING2 protein 1-280 amino acids purified from E. coli.
Description	Inhibitor of Growth Family, Member 1 Human Recombinant

Target information (Q9UK53)

Gene Name

About Gene Symbol