

ISCU Human

Iron-sulfur Cluster Scaffold Homolog Human Recombinant

RCP1132

Product Overview

Name	ISCU Human
Accession(Primary)	Q9H1K1
Immunogen	Anti-human ISG15 mAb is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human ISG15 amino acids 1-157 purified from E. coli.
Description	Iron-sulfur Cluster Scaffold Homolog Human Recombinant

Target information (Q9H1K1)

Gene Name**About Gene Symbol**