

ISG15 Human

ISG15 Human Recombinant

RCP1133

Product Overview

Name	ISG15 Human
Accession(Primary)	P05161
Description	ISG15 Human Recombinant

Target information (P05161)

Gene Name

About Gene Symbol

ISG15 (ISG15 Ubiquitin Like Modifier) is a cytokine and immune signaling protein involved in JAK-STAT/cytokine signaling and ubiquitin-proteasome pathway. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.