

KCTD11 Human

Potassium Channel Tetramerisation Domain Containing 11 Human Recombinant

RCP1153

Product Overview

Name	KCTD11 Human
Accession(Primary)	Q693B1
Immunogen	Anti-human KCTD15 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human KCTD15 amino acids 1-234 purified from E. coli.
Description	Potassium Channel Tetramerisation Domain Containing 11 Human Recombinant

Target information (Q693B1)

Gene Name

About Gene Symbol