

KCTD15 Human

Potassium Channel Tetramerisation Domain Containing 15 Human Recombinant

RCP1154

Product Overview

Name	KCTD15 Human
Accession(Primary)	Q96SI1
Immunogen	Anti-human KCTD15 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human KCTD15 protein 1-159 amino acids purified from E. coli.
Description	Potassium Channel Tetramerisation Domain Containing 15 Human Recombinant

Target information (Q96SI1)

Gene Name

About Gene Symbol