

KIR2DL1 Human

Killer Cell Immunoglobulin-Like Receptor, 2 Domains Long Cytoplasmic Tail 1 Human

RCP1161

Product Overview

Name	KIR2DL1 Human
Accession(Primary)	P43626
Immunogen	Anti-human KIR2DL3 mAb, is derived from hybridization of mouse P3-x63-Ag8.653 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human KIR2DL3 amino acids 19-161 purified from E. coli.
Description	Killer Cell Immunoglobulin-Like Receptor, 2 Domains Long Cytoplasmic Tail 1 Human Recombinant

Target information (P43626)

Gene Name

About Gene Symbol