

KIR2DL4 Human

Killer Cell Immunoglobulin-Like Receptor, 2 Domains Long Cytoplasmic Tail 4 Human

RCP1163

Product Overview

Name	KIR2DL4 Human
Accession(Primary)	Q99706
Immunogen	Anti-human KIR2DL5A mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human KIR2DL5A amino acids 22-238 purified from E. coli.
Description	Killer Cell Immunoglobulin-Like Receptor, 2 Domains Long Cytoplasmic Tail 4 Human Recombinant

Target information (Q99706)

Gene Name

About Gene Symbol