

KIR2DL5A Human

Killer Cell Immunoglobulin-Like Receptor, 2 Domains Long Cytoplasmic Tail 5A Human

RCP1164

Product Overview

Name	KIR2DL5A Human
Accession(Primary)	Q8N109
Immunogen	Anti-human KIR2DS4 mAb is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human KIR2DS4 purified from E. coli.
Description	Killer Cell Immunoglobulin-Like Receptor, 2 Domains Long Cytoplasmic Tail 5A Human Recombinant

Target information (Q8N109)

Gene Name

About Gene Symbol