

KLRC3 Human

Killer Cell Lectin-Like Receptor Subfamily C, Member 3 Human Recombinant

RCP1181

Product Overview

Name	KLRC3 Human
Accession(Primary)	Q07444
Immunogen	Anti-human KLRD1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human KLRD1 protein 32-279 amino acids purified from E. coli.
Description	Killer Cell Lectin-Like Receptor Subfamily C, Member 3 Human Recombinant

Target information (Q07444)

Gene Name

About Gene Symbol