

KRT20 Human, His

Cytokeratin 20 Human Recombinant, His Tag

RCP1200

Product Overview

Name	KRT20 Human, His
Accession(Primary)	P35900
Immunogen	Anti-human KRT23 mAb, is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice immunized with recombinant human KRT23 amino acids 271-422 purified from E. coli.
Description	Cytokeratin 20 Human Recombinant, His Tag

Target information (P35900)

Gene Name

About Gene Symbol