

L1CAM Human

L1 Cell Adhesion Molecule Human Recombinant

RCP1209

Product Overview

Name	L1CAM Human
Accession(Primary)	P32004
Description	L1 Cell Adhesion Molecule Human Recombinant

Target information (P32004)

Gene Name

About Gene Symbol

L1CAM (L1 Cell Adhesion Molecule) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.