

LAMP2 Human

Lysosomal-Associated Membrane Protein 2 Human Recombinant

RCP1219

Product Overview

Name	LAMP2 Human
Accession(Primary)	P13473
Description	Lysosomal-Associated Membrane Protein 2 Human Recombinant

Target information (P13473)

Gene Name

About Gene Symbol

LAMP2 (Lysosomal Associated Membrane Protein 2) is a cell surface marker involved in adaptive immunity and vesicular trafficking. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.