

LILRA3 Human

Leukocyte Immunoglobulin Like Receptor A3 Human Recombinant

RCP1240

Product Overview

Name	LILRA3 Human
Accession(Primary)	Q8N6C8
Description	Leukocyte Immunoglobulin Like Receptor A3 Human Recombinant

Target information (Q8N6C8)

Gene Name

About Gene Symbol

Background

Leukocyte Immunoglobulin Like Receptor A3 (LILRA3) acts as a soluble receptor for class I MHC antigens. LILRA3 binds classical and non-classical HLA class I molecules and exhibits high affinity for monocyte surfaces, effectively suppressing LPS-induced TNF-alpha production. Unlike many of its family, LILRA3 has no transmembrane domain.