

LRP4 Human

LDL Receptor Related Protein 4 Human Recombinant

RCP1264

Product Overview

Name	LRP4 Human
Accession(Primary)	O75096
Immunogen	Anti-human LRPAP1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human LRPAP1 protein 35-357 amino acids purified from E. coli.
Description	LDL Receptor Related Protein 4 Human Recombinant

Target information (O75096)

Gene Name

About Gene Symbol