

LSM12 Human

LSM12 Homolog, U6 Small Nuclear RNA Associated Human Recombinant

RCP1270

Product Overview

Name	LSM12 Human
Accession(Primary)	Q3MHD2
Immunogen	Anti-human LSM2 mAb, is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice immunized with recombinant human LSM2 amino acids 1-95 purified from E. coli.
Description	LSM12 Homolog, U6 Small Nuclear RNA Associated Human Recombinant

Target information (Q3MHD2)

Gene Name

About Gene Symbol