

LSM4 Human

LSM4 Homolog, U6 Small Nuclear RNA Associated Human Recombinant

RCP1273

Product Overview

Name	LSM4 Human
Accession(Primary)	Q9Y4Z0
Immunogen	Anti-human LSM5 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human LSM5 protein 1-91 amino acids purified from E. coli.
Description	LSM4 Homolog, U6 Small Nuclear RNA Associated Human Recombinant

Target information (Q9Y4Z0)

Gene Name

About Gene Symbol