

LSM5 Human

LSM5 Homolog, U6 Small Nuclear RNA Associated Human Recombinant

RCP1274

Product Overview

Name	LSM5 Human
Accession(Primary)	Q9Y4Y9
Immunogen	Anti-human LTA mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human LTA amino acids 35-205 purified from E. coli.
Description	LSM5 Homolog, U6 Small Nuclear RNA Associated Human Recombinant

Target information (Q9Y4Y9)

Gene Name

About Gene Symbol