

## MFAP3 Human

*Microfibrillar-associated Protein 3 Human Recombinant*

RCP1361

### Product Overview

|                           |                                                                                                                                                                                                                               |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>               | MFAP3 Human                                                                                                                                                                                                                   |
| <b>Accession(Primary)</b> | P55082                                                                                                                                                                                                                        |
| <b>Immunogen</b>          | Anti-human MFAP4 mAb, clone PAT12D11AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human MFAP4 protein 22-255 amino acids purified from E. coli. |
| <b>Description</b>        | Microfibrillar-associated Protein 3 Human Recombinant                                                                                                                                                                         |

### Target information (P55082)

**Gene Name**

**About Gene Symbol**