

MRPS23 Human

Mitochondrial Ribosomal Protein S23 Human Recombinant

RCP1401

Product Overview

Name	MRPS23 Human
Accession(Primary)	Q9Y3D9
Immunogen	Anti-human MRPS25 mAb, clone PAT38E7A, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human MRPS25 protein 1-173 amino acids purified from E. coli.
Description	Mitochondrial Ribosomal Protein S23 Human Recombinant

Target information (Q9Y3D9)

Gene Name

About Gene Symbol