

MRPS28 Human

Mitochondrial Ribosomal Protein S28 Human Recombinant

RCP1403

Product Overview

Name	MRPS28 Human
Accession(Primary)	Q9Y2Q9
Immunogen	Anti-human MRRF mAb, clone PAT7D10A, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human MRRF protein 56-262 amino acids purified from E. coli.
Description	Mitochondrial Ribosomal Protein S28 Human Recombinant

Target information (Q9Y2Q9)

Gene Name

About Gene Symbol