

MYL2 Human

Myosin Light Chain 2 Human Recombinant

RCP1429

Product Overview

Name	MYL2 Human
Accession(Primary)	P10916
Immunogen	Anti-human MYL4 mAb, clone PAT4E8AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human MYL4 protein 1-197 amino acids purified from E. coli.
Description	Myosin Light Chain 2 Human Recombinant

Target information (P10916)

Gene Name**About Gene Symbol**