

NKIRAS2 Human

NFKB Inhibitor Interacting Ras-Like 2 Human Recombinant

RCP1499

Product Overview

Name	NKIRAS2 Human
Accession(Primary)	Q9NYR9
Immunogen	Anti-human NKp30 mAb, is derived from hybridization of mouse P3-X63-Ag8.653 myeloma cells with spleen cells from BALB/c immunized with recombinant human NKp30 amino acids 19-138 purified from E. coli.
Description	NFKB Inhibitor Interacting Ras-Like 2 Human Recombinant

Target information (Q9NYR9)

Gene Name

About Gene Symbol