

PABPN1 Human

Poly(A) Binding Protein, Nuclear 1 Human Recombinant

RCP1581

Product Overview

Name	PABPN1 Human
Accession(Primary)	Q86U42
Description	Poly(A) Binding Protein, Nuclear 1 Human Recombinant

Target information (Q86U42)

Gene Name

About Gene Symbol

PABPN1 (Poly(A) Binding Protein Nuclear 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.