

PARK7 Mouse

Parkinson Disease Protein 7 Mouse Recombinant

RCP1588

Product Overview

Name	PARK7 Mouse
Accession(Primary)	Q99LX0
Immunogen	Anti-human PARK7 mAb, clone P1E12AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PARK7 protein 1-189 amino acids purified from E. coli
Description	Parkinson Disease Protein 7 Mouse Recombinant

Target information (Q99LX0)

Gene Name

About Gene Symbol