

PCMTD1 Human

Protein-L-Isoaspartate O-Methyltransferase 1 Human Recombinant

RCP1596

Product Overview

Name	PCMTD1 Human
Accession(Primary)	Q96MG8
Immunogen	Anti-human PCNA mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human PCNA amino acids 1-261 purified from E. coli.
Description	Protein-L-Isoaspartate O-Methyltransferase 1 Human Recombinant

Target information (Q96MG8)

Gene Name

About Gene Symbol